

GRANULES & CUBES

Granules & cubes, cancellous

Origin: Human
 Tissue: Cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Indication: Bone void filler
 Rehydration: Min. 10 minutes

Granules, cortico-cancellous

Origin: Human
 Tissue: Cortico-cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Indication: Bone void filler
 Rehydration: Min. 10 minutes

GRAIN SIZES OF GRANULES

The granule sizes are achieved by sieving. The different perforation of the sieves leads to the sizes listed below. Depending on the direction of fall, particles in one dimension can be slightly larger than specified.

ORDER INFORMATION*



DESCRIPTION	GRAIN SIZE [mm]	ITEM NUMBER	VOLUME [cc]
Cancellous granules	2 - 5	ALO319	5
		ALO315	15
		ALO309	30
		ALO317	45
	5 - 8	ALO326	5
		ALO316	15
		ALO310	30
		ALO331	45
	2 - 8	ALO305	5
		ALO306	15
ALO307		30	
ALO308		45	
> 8	ALO300	15	
	ALO301	30	
Cancellous granules - Spierings	5 - 10	ALO350	30
	2 - 8	ALO351	10
ALO352		15	
ALO353		30	
Cortico-cancellous granules	2 - 8	ALO340	15
		ALO341	30
Cancellous cubes*	5x5x5	ALO325	10
		ALO314	20

* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

