

DEMINERALIZED BONE MATRIX - DBM

DBM Putty

Origin: Deceased donor, human
Tissue: 100% allogenic bone
Processing: Demineralization
Sterilisation: Gamma irradiation
Indication: Bone void filler

DBM Fibers

Origin: Deceased donor, human
Tissue: 100% allogenic bone, cortical
Processing: Demineralization/Lyophilisation
Sterilisation: Gamma irradiation
Indication: Bone void filler

DBM Putty

DESCRIPTION	ITEM NUMBER	VOLUME [cc]
Maxxeus DBM Putty	CTS2666 - 40	1,0
	CTS2667 - 40	2,5
	CTS2668 - 40	5,0
	CTS2669 - 40	10,0



Maxxeus DBM Putty

DBM Fibers

DESCRIPTION	ITEM NUMBER	VOLUME [cc]
Maxxeus DBM Fibers	CTS3594 - 12	1,0
	CTS3595 - 12	2,5
	CTS3596 - 12	5,0
	CTS3597 - 12	10,0
	CTS3598 - 12	15,0

Maxxeus DBM Fibers

