


FEMORAL HEAD

Femoral head, halved

Origin: Human
 Tissue: Cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Indication: Bone void filler
 Rehydration: Min. 10 minutes

Femoral head, bisected


Origin: Human
 Tissue: Cortico-cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Indication: Bone void filler
 Rehydration: Min. 10 minutes



ORDER INFORMATION*




| DESCRIPTION | ITEM NUMBER | D [mm] | H [mm] |
|----------------------|-------------|--------|--------|
| Femoral head, halved | ALO441 | < 45 | 20 |
| | ALO444 | > 45 | 20 |



ORDER INFORMATION*



| DESCRIPTION | ITEM NUMBER | SIZE |
|------------------------|-------------|-------|
| Femoral head, bisected | ALO446 | short |
| | ALO447 | long |



* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.